
Dental Cement Classification and Properties

Classification	Color Key	Product Name		Adhesive Required
Glass lonomer		Ketac™ Cem Cement		No
Resin Modified Glass lonomer		RelyX™ Luting Cement	A STATE OF THE PARTY OF THE PAR	No
		RelyX [™] Luting Plus Cement	Part of the second	No
Resin		RelyX [™] Unicem Self-Adhesive Cement		No
		RelyX [™] ARC Cement	\$ my	Yes
		RelyX™ Veneer Cement		Yes

Properties









^{*} Costs are based on estimates per application and are subject to change. Information is presented as a generalization.



Cement Products – Recommended Indications*

Restoration Type	RelyX™ Unicem Self-Adhesive Universal Resin Cement	RelyX™ Luting RelyX™ Luting Plus Cements	RelyX™ Veneer Cement	RelyX™ ARC Adhesive Resin Cement	Ketac™ Cem Radiopaque Ketac™ Cem Aplicap™/ Maxicap™ Permanent Glass Ionomer Luting Cements
Metal/PFMs	-	_		_	
Porcelain/Ceramic					
Veneers					
Zirconia					
CAD/CAM					
Posts					
Pediatric Crowns • Metal • Composite		<u>_</u>			
Maryland Bridges				_	

^{*}These are only recommended indications. For a complete list, refer to the respective product Instructions for Use.

Delivery Method





Clicker[™] Dispenser



Powder/Liquid



Ordering Information



Product Information Item # RelyX[™] Unicem Self-Adhesive Universal Resin Cement in the Clicker™ Dispenser **Single Packs** Kit Includes: 1 Clicker Dispenser**, 1 Mix Pad,

Instructions, Technique Card

56838 **Shade Translucent** 56837 **Shade A2 Universal** 56836 Shade A3 Opaque

Additional RelyX Unicem cement packages and accessories are available in Aplicap™ and Maxicap™ capsule delivery systems. Ask your certified 3M ESPE distributor for more information.



Item # **Product Information**

RelyX[™] Luting Plus Cement Trial Kit

Kit Includes: 1 RelyX™ Luting Plus Cement in Clicker™ Dispenser**, Illustrated Technique Guide and Instructions for Use

RelyX[™] Luting Plus Cement Refills

Kit Includes: 2 RelyX $^{\text{\tiny{TM}}}$ Luting Plus Cement in Clicker $^{\text{\tiny{TM}}}$ Dispensers ** and Instructions for Use

^{**}Each Clicker dispenser contains 11.0 grams of material, dispensed in 80 clicks (approximately 40 applications)



Customer Care Center: 1-800-634-2249 www.3MESPE.com/relyXcements

Dental Products

3M Center Building 275-2SE-03 St. Paul, MN 55144-1000 U.S.A. 1-800-634-2249

3M Canada Post Office Box 5757 London, Ontario N6A 4T1 Canada 1-888-363-3685