


KATANA™ Zirconia



New High Esthetic Potential for Zirconia in Dental Restorations *

NEW


UTML Ultra Translucent Multi Layered

STML Super Translucent Multi Layered

ML Multi Layered **HT** High Translucent

The translucency is 1.4 times more than the conventional zirconia products.



The translucency of zirconia raw materials was increased 1.4 times by our proprietary research and development. Now you can achieve natural life like translucency with KATANA™ monolithic zirconia materials.



Date source: Kuraray Noritake Dental Inc. The numerical value varies according to a condition

You now have a wide range of clinical case applications for monolithic zirconia with the new additional KATANA™ series.

Selecting the right translucency and mechanical properties, you can fabricate many options, from the esthetic anterior restorations to the posterior bridge work.*1



Date source: Kuraray Noritake Dental Inc. The numerical value varies according to a condition

*1 UTML, STML, ML : The minimum thickness of zirconia which is necessary in each series is different. Read carefully the instructions prior to using the material.
*2 HT is recommended for the frame-work if you overlay with layered porcelains.